iCD. GmbH
Augustinusstr. 9d
50226 Frechen / Germany
Phone +49 2234 96634-0
Fax +49 2234 96634-90
Email: info@icd.eu
Web: www.icd.eu

LABS/Q® Actions and measures

LABS/Q® Actions and measures

In a laboratory many activities occur which have to be carried out in a defined workflow according to specifications, SOPs, work or process instructions. There may be processes like OOS (out of specification), OOT (out of trend) or the handling of HACCP tests (Hazard Analysis and Critical Control Point). It could also be a DIN ISO 9000 audit with deviations occurring and ensuing actions and measures.


LABS/Q®-version 3.6 has implemented new, optional modules and functions to determine workflow actions and measures as well as carrying them out and evaluating data. Subsequently you will find a fictitious OOS procedure which shows how the newly implanted functions of **LABS/Q®** support laboratory personnel.

Planning a workflow to process an OOS

Company procedures are determined by responsibilities and company structures as well as existing SOPs. A special module is available to represent a specific workflow in **LABS/Q®**. This module determines how an OOS action and its measures are processed.

These are the stipulations:

- **Changes of status**
Chang-over of status to OOS as well as to its measures can be defined individually. You can determine which switch from status to status is possible.
- **Following actions**
After causal research is finished two repetition operations have to be done. Or a complete method repetition is necessary. Or a new sample has to be processed. Maybe a training scheme or a control of instrument or control type is necessary.
- **Information**
A distribution list for the OOS action and its possible measures is determined. It contains information who or which group of people will receive an automatic report following a specific action.
- **Responsibilities**
It can be defined who is responsible for which OOS action or following measures and how these persons can be automatically determined by the system.
- **Date targets**
A maximum period to process an action is determined here. Based on these specifications a deadline monitoring can be effected manually or automatically.
- **Check lists for individual actions**
Check lists and anticipated values can be determined here. Further procedures are initiated according to the results of the check lists.

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
Opening an OOS action

In case of a specification violation the deviation must be dealt with as an OOS. Depending on the particular processing (e.g. sample sequence during chromatography) 1 up to n samples may be affected by this action.

LABS/Q® offers the option to include an OOS action during all results recording functions. All affected properties, even cross-sample ones, can be selected and assigned to the OOS. Depending on the configuration of **LABS/Q®** a submission release is only possible when the complete OOS is processed starting with the opening of an OOS action.


Example: Opening an OOS action within a result entry

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The OOS with a first measure is processed based on the specification “type of action” and “level of work” as well as planning information for the OOS and its actions.

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Processing an OOS and its measures

In case an OOS is pending its execution is often regulated via an SOP. Often an internal laboratory causal research is conducted first. This makes use of check lists as tools. Depending on the results authorised persons determine whether a repeated operation or a complete methodical repetition must be done or if even a new sample has to be analysed. After finishing the analysis the OOS is either closed or further actions are taken, e.g. a training session for employees. Maybe the method has to be revised.

In **LABS/Q®** the first measure (e.g. laboratory causal research) can be executed. When the measure is finished consequent measures can be added. You can determine for each measure whether the submission


Example: OOS action with 1 action and a check list

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has to be supplemented by further operations or by additional samples. In case of the latter the submission is automatically supplemented by these specifications and the adequate status.

Example: Inserting further actions to the OOS

This screenshot shows the 'Maßnahme zum OOS einfügen' dialog box. It allows adding a new measure with fields for 'Vorlagen ID' (MASSNAHME), 'Aktionstyp' (Out of Spec.), and 'Bearbeitungsstufe' (Stufe 1). The 'Aktions-Informationen' section on the right shows the current OOS details, including the 'Bemerkung' and 'SOP-Dokument' fields. A table at the bottom shows the history of the OOS, listing creation and modification dates and users.

Evaluation of OOS/OOT actions

In **LABS/Q**® appropriate evaluation functions are at your disposal. These evaluations can be designed individually.

Example: Survey over all OOS's/OOTs during a period of time

HTML-View and Print - [c:\icd\labsq350\labsq350\00158908.htm]

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icd.GmbH